
Bambach Saddle Seat®

Measurement Guide





Height From The Seat To The Floor


Short Cylinder: lowest setting 20" (51cm) / highest setting 25" (63cm)
Medium Cylinder: lowest setting 24" (60cm) / highest setting 30" (77cm)
Tall Cylinder: lowest setting 26" (66cm) / highest setting 36" (91cm)

Seat Size Guidance

Classic: Men up to 36" waist and women size 12 -16 waist.

Classic Plus: Men size 38" waist and over and woman size 18 waist and over.

Classic Small: Woman up to size 8-10 waist and children.




Bambach Saddle Seat®

How To Choose Your Cylinder Height










User Height: 5' 2" and under








User Height: 5′ 3″ to 6′





User Height: 6' 1" and over



^{*} When a foot ring is being used, such as with a dental assistant or drafting stool, a tall cylinder is typically recommended for most user heights